


KIMAX[®] Heavy-Duty Glassware ...

tough enough

for your front line.



KIMAX[®] borosilicate glassware for maximum mechanical strength and thermal shock resistance in your hard-hitting laboratory conditions.



888-546-2531
www.kimble-chase.com
info@kimble-chase.com

Try these flasks
and you'll become
a fan for life!



92812N

92812G

See the web for
flasks with other
stopper options:
www.kimble-chase.com

PTFE



PE



Choose
RAY-SORB®
processed
flasks to
protect your
light-sensitive
contents.



92822N

92822G

KIMAX® Heavy-Duty Wide-Mouth Volumetric Flasks

- Superior strength, durability and safety from heavy, uniform walls
- Wide-mouth facilitates up to 50mL pipet access
- Quick identification with large, permanent, volume markings
- Flat interior bottoms for stirring with standard stir bars
- Accuracy - graduated to Class A volumetric tolerances
- Choose RAY-SORB® for light-sensitive compounds such as Vitamin A

KIMAX® Volumetric Flasks with Glass ⌘ Stoppers

Catalog Number	Cap. (mL)	±Tol. (mL)	⌘ Stopper	Approx. Height with Stopper (mm)	Case Qty.	Case Price
92812G-5	5	0.08	13	97	6	144.59
92812G-10	10	0.08	13	100	6	155.87
92812G-20	20	0.08	13	130	6	170.46
92812G-25	25	0.08	13	130	6	167.66
92812G-50	50	0.08	13	150	6	172.10
92812G-100	100	0.10	16	180	6	209.58
92812G-200	200	0.20	19	223	6	272.45
92812G-250	250	0.20	19	248	6	272.45
92812G-500	500	0.20	19	283	6	326.94
92812G-1000	1000	0.30	22	335	1	102.12
92812G-2000	2000	0.50	27	395	1	149.52

KIMAX® Volumetric Flasks without Stoppers

Catalog Number	Cap. (mL)	±Tol. (mL)	⌘ Stopper Size	Approx. Height without Stopper (mm)	Case Qty.	Case Price
92812N-5	5	0.08	13	77	6	135.77
92812N-10	10	0.08	13	87	6	147.06
92812N-20	20	0.08	13	110	6	161.72
92812N-25	25	0.08	13	110	6	158.93
92812N-50	50	0.08	13	130	6	163.28
92812N-100	100	0.10	16	160	6	200.84
92812N-200	200	0.20	19	200	6	263.72
92812N-250	250	0.20	19	225	6	263.72
92812N-500	500	0.20	19	260	6	318.21
92812N-1000	1000	0.30	22	310	1	101.63
92812N-2000	2000	0.50	27	370	1	148.05

KIMAX® RAY-SORB® Volumetric Flasks with Glass ⌘ Stoppers

Catalog Number	Cap. (mL)	±Tol. (mL)	⌘ Stopper	Approx. Height with Stopper (mm)	Case Qty.	Case Price
92822G-5	5	0.08	13	97	6	179.55
92822G-10	10	0.08	13	100	6	192.81
92822G-20	20	0.08	13	130	6	198.41
92822G-25	25	0.08	13	130	6	198.41
92822G-50	50	0.08	13	150	6	214.47
92822G-100	100	0.10	16	180	6	291.66
92822G-200	200	0.20	19	223	6	321.63
92822G-250	250	0.20	19	248	6	338.55
92822G-500	500	0.20	19	283	6	398.51
92822G-1000	1000	0.30	22	335	1	124.76
92822G-2000	2000	0.50	27	395	1	180.76

KIMAX® RAY-SORB® Volumetric Flasks without Stoppers

Catalog Number	Cap. (mL)	±Tol. (mL)	⌘ Stopper Size	Approx. Height without Stopper (mm)	Case Qty.	Case Price
92822N-5	5	0.08	13	77	6	170.81
92822N-10	10	0.08	13	80	6	184.08
92822N-20	20	0.08	13	110	6	189.68
92822N-25	25	0.08	13	110	6	189.68
92822N-50	50	0.08	13	130	6	205.74
92822N-100	100	0.10	16	160	6	282.93
92822N-200	200	0.20	19	200	6	312.82
92822N-250	250	0.20	19	225	6	326.59
92822N-500	500	0.20	19	260	6	385.98
92822N-1000	1000	0.30	22	310	1	123.30
92822N-2000	2000	0.50	27	370	1	179.31

KIMAX®




When safety is a concern, you can gain some extra yardage with heavy-duty glassware from Kimble Chase.

KIMAX® Heavy-Duty Low Form Beakers with Double Capacity Scales

- Made from heavy-wall tubing for uniform wall thickness
- Uniformity of construction allows for use on hot plates
- Evenly tooled top rim and spout to facilitate pouring
- Durable matte finish marking area and white graduations
- Designed form ASTM Specification E960, Type II requirements

Thick, uniform walls for increased resistance to mechanical shock.



14005

Catalog Number	Capacity (mL)	Graduated Interval (mL)	Subdivision (mL)	Approx. OD x Height (mm)	Shelf-Pack Qty.	Case Qty.	Case Price
14005-250	250	25 to 200	25	67 x 92	12	48	311.51
14005-400	400	25 to 325	25	77 x 114	12	48	406.34
14005-600	600	50 to 500	50	87 x 122	6	36	355.18
14005-1000	1000	100 to 1000	50	107 x 152	6	24	451.75
14005-2000	2000	200 to 1800	200	130 x 190	4	8	307.77
14005-4000	4000	500 to 3500	500	161 x 252	1	4	300.54

KIMAX® Heavy-Duty Solution Bottle (Carboy)

- Neck is tooled for a #12 stopper
- 5 gallon size is designed from Federal Specification DD-B-597

Catalog Number	Approx. Capacity Liters - Gallons	Approx. Dia. (mm)	Approx. Height (mm)	Case Qty.	Case Price
14950-25	9.5 - 2.5	222	392	1	160.34
14950-35	13.2 - 3.5	257	448	1	209.27
14950-500	19.0 - 5.0	294	502	1	314.45
14950-120	45.5 - 12.0	410	584	1	705.89



14950

Carboy Clamp and Silicone Stopper Assemblies for 14950 Solution Bottles

- Two or Three-port, platinum cured stopper assembly
- Sanitary KYNAR® clamp for sterile transfer with no extractables
- Ports connect to 1/8", 1/4" (4-port stopper has two 1/4" tubes) and 3/8" with 2' silicone tubing at both ends



Accessory Clamp and Stopper Assembly

Catalog Number	Description	Tubing Connections	Case Qty.	Case Price
14950C-12	Clamp	-	1	44.64
14950S-321	Stopper, 3-port	1/8"(1) 1/4"(1) 3/8"(1)	1	223.99
14950S-2321	Stopper, 4-port	1/8"(1) 1/4"(2) 3/8"(1)	1	254.65

KYNAR is a registered trademark of ARKEMA for polyvinylidene fluoride.





Try KimCote plastic-coated glassware for an added measure of safety!

- Helps prevent and contain spills in the event of labware breakage
- Ultra-clear for improved sample viewing
- Non-slip surface — wet or dry
- Autoclavable and chemical resistant
- Durable and odor-free



For more information: 888-546-2531 • www.kimble-chase.com • info@kimble-chase.com

Make the right play with KIMAX® heavy-duty glassware.



See the web for flasks with other hose connection options:

www.kimble-chase.com



27060

KIMAX® Heavy-Wall Filter Flasks with Graduations

- Designed to hold 1 atmosphere of vacuum or 29 inches of mercury
- Sizes 25 - 500 mL have hose connections for 5/16-inch I.D. flexible tubing
- Sizes 1000 - 4000 mL have hose connections for 3/8-inch I.D. flexible tubing
- Designed from ASTM Specification E1406, Type II, Class I requirements

Catalog Number	Cap. (mL)	Graduated Interval (mL)	Subdiv. (mL)	Rubber Stopper Size	Approx. Widest Dia. (mm)	Approx. Height (mm)	Shelf-Pack Qty.	Case Qty.	Case Price
27060-25	25	5 to 25	5	3	40	74	6	18	396.94
27060-50	50	20 to 50	10	4	48	81	6	18	396.94
27060-125	125	50 to 125	25	5	69	111	6	18	410.71
27060-250	250	50 to 250	25	6	83	163	6	18	305.17
27060-500	500	150 to 500	50	7	104	190	6	18	529.32
27060-1000	1000	300 to 1000	50	8	135	238	6	12	556.92
27060-2000	2000	600 to 2000	100	9	168	300	—	1	92.18
27060-4000	4000	1000 to 4000	250	12	208	380	—	1	135.39



45212



45196


See our online catalog for a variety of centrifuge tubes.

KIMAX® Heavy-Duty Screw Thread Centrifuge Tubes with Caps

- Made from heavier wall tubing to withstand higher centrifuge speeds
- Choose from tapered or round bottoms
- Phenolic caps with cemented-in rubber liners (45066B) supplied but not attached

Catalog Number	Description	Capacity (mL)	O.D. x Length (mm)	GPI Thread Finish	Case Qty.	Case Price
45196-40	tapered bottom	40	29 x 140	24-410	12	290.02
45212-35	round bottom	35	29 x 100	24-410	12	142.65
45212-50	round bottom	50	29 x 123	24-410	12	205.43

www.kimble-chase.com



We're here to help.

If you don't find what you need, please call for assistance. Our customer service group will be glad to discuss your special requirements.

To order, contact your local distributor or the office listed below.
Customer Service Hours: 8 am - 5 pm EST

Kimble Chase

1022 Spruce Street
PO Box 1502
Vineland NJ 08362-1502

Toll-Free: 888-546-2531
www.kimble-chase.com
info@kimble-chase.com



Printed in USA. Trademarks are the property of their respective owners. Specifications and prices subject to change without notice.
Prices shown are US dollars.

RJ5Md5Me0611KC1104